Product: **ProTaper® Universal Retreatment Kit**
REF A1410 – A1412 – A1415

Ref doc.: **F19 02 89.X - F19 02 86.X/ 04 / 2006 updated 04/2009**

		PAGE
GB	DIRECTIONS FOR USE (ProTaper® Universal Retreatment Kit)	2
F	PROTOCOLE D'UTILISATION (ProTaper® Universal Retreatment Kit)	7
D	GEBRAUCHSANLEITUNG (ProTaper® Universal Retreatment Kit)	12
I	MODO DI USO (ProTaper® Universal Retreatment Kit)	17
E	MODO OE UTILIZACION (ProTaper® Universal Retreatment Kit)	22
NL	GEBRUIKSVOORSCHRIFT (ProTaper® Universal Retreatment Kit)	27
DK	BRUGSANVISNING (ProTaper® Universal Retreatment Kit)	32
FIN	KÄYTTÖOHJEET (ProTaper® Universal Retreatment Kit)	37
S	BRUKSANVISNING (ProTaper® Universal Retreatment Kit)	42
P	PROTOCOLO DE UTILIZAÇÃO (ProTaper® Universal Retreatment Kit)	47
GR	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ ΓΙΑ (ProTaper® Universal Retreatment Kit)	52
POL	INSTRUKCJA OBSŁUGI (ProTaper® Universal Retreatment Kit)	57

For Dental Use Only

DIRECTIONS FOR USE (ProTaper® Universal Retreatment Kit)

PROTAPER® INSTRUMENTS FOR ROOT CANAL RETREATMENT

- ProTaper® D1: file to unfill the coronal one third
- ProTaper® D2: file to unfill the middle one third
- ProTaper® D3: file to unfill the apical one third

COMPOSITION: the cutting part of these instruments is made of a nickel-titanium alloy.

1) INDICATIONS FOR USE

These instruments are to be used only in a clinical or hospital environment, by qualified users.

Application field: canal unfilling for endodontic retreatment.

2) CONTRAINDICATIONS

- Do not use to remove resin type pastes.
- Do not use the D1 file around a canal curvature.

3) WARNINGS

- ProTaper® retreatment files are used to unfill root canals filled with gutta percha points, Thermafil obturators or, eugenol based soluble paste. They can't be used to unfill resine type paste.

4) PRECAUTIONS

- Multiple use disinfection and resterilization cycles may lead to increased risk of file separation.
- These instruments should not be immersed in a sodium hypochlorite solution.
- Instrument decontamination. strictly follow decontamination instructions from the manufacturer.
- Use in a constant rotation speed of:
 - 500 rpm for removing gutta percha or Thermafil obturators
 - 250 – 300 rpm for removing zinc oxide – eugenol based soluble paste.
- For optimal usage, torque control devices are recommended.
- Use the D1 file for the coronal one third of the canal, then the D2 file for the middle one and at last the D3 file for the apical one third.

5) ADVERSE REACTIONS

In the present technical state, no adverse reaction has been reported so far.

6) STEP BY STEP INSTRUCTIONS FOR PROTAPER FILES

Before removing gutta percha, carrier-based obturators or paste from a root canal:

- Carefully observe 3 different, horizontally angulated radiographic images.
- Visualize the density of obturation material relative to the width, length, and curvature of the canal.
- Access the pulp chamber and note the circumferential dimensions of the obturation material at the orifice(s).
- Select the best removal technique after radiographic and clinical assessment.
- Without cutting dentin, remove obturation material in a progressive crown-down manner.

Gutta Percha / Carrier-Based Obturator Removal

- When the rotary removal method is utilized, select the lowest speed (500 RPM) that will effectively engage and remove obturation material from the canal.
- Without engaging dentin, gently press the spinning ProTaper® D1 file into the gutta percha to create friction, generate a heat wave, and extract material out of the canal. Never engage ProTaper® D1 around a canal curvature.*
- Remove the D1 file frequently, inspect the blades for obturation material and clean the debris from the flutes.
- Continue with the D1 file, or the ProTaper® Retreatment file that fits passively between the dentinal walls, until gutta percha is removed from the coronal one-third of the canal.
- Select the ProTaper® D2 file and, using one or more passes, extract obturation material from the middle one-third of the canal. Use a brushing outstroke motion to remove material from the canal walls.
- When appropriate, choose the ProTaper® D3 file and lightly press into the more deeply positioned material and auger obturation material out of the apical one-third of the canal.
- Continue with the D3 file as long as the flutes of the instrument, upon removal, are loaded with obturation material.
- When the obturation material is short of the canal terminus, use small sized hand files in the presence of a viscous chelator to negotiate and secure the rest of the canal.
- After assessing the glide path, select either manual or rotary NiTi ProTaper® files to shape and finish the canal.

* In the instance of carrier removal, select the appropriately tapered ProTaper® Retreatment file that can be carried sufficiently deep into the canal and lateral to the carrier. A long engagement zone will more effectively auger the entire length of the carrier out of the canal.

Soluble Zinc oxide – Eugenol based Paste Removal

- When the rotary removal method is utilized, select the lowest speed (250-300 RPM) that will effectively engage and remove obturation material from the canal.
- Flood the pulp chamber with the appropriate solvent (DMS IV from Dentsply Maillefer) and probe the canal orifice with an explorer to check if the paste has been effectively softened.
- Without engaging dentin, gently press the spinning ProTaper® D1 file into the material and use a short pecking motion to extract material out of the canal. Never engage ProTaper® D1 around a canal curvature.
- Remove the D1 file frequently, inspect the blades for obturation material and clean the debris from the flutes.
- Continue with the D1 file, or the ProTaper® Retreatment file that fits passively between the dentinal walls, until paste is removed from the coronal one-third of the canal.
- Select the ProTaper® D2 file and repeat the same pecking action to extract obturation material from the middle one-third of the canal. Use a brushing outstroke motion to remove material from the canal walls.
- When appropriate, choose the ProTaper® D3 file and, in the same way auger the more deeply positioned paste material out of the apical one-third of the canal.
- Continue with the D3 file as long as the flutes of the instrument, upon removal, are loaded with obturation material.
- When the obturation material is short of the canal terminus, use small sized hand files in the presence of a viscous chelator to negotiate and secure the rest of the canal.
- After assessing the glide path, select either manual or rotary NiTi ProTaper® files to shape and finish the canal.

7) DISINFECTION, CLEANING AND STERILIZATION

Reprocessing procedure for dental instruments and implantable radicular devices

Foreword

For hygiene and sanitary safety purposes, all instruments not marked "sterile" must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well as the subsequent ones.

Area of application

Disinfection and sterilisation before first usage and reprocessing procedures concerning:

A1. Instruments

Cutting instruments, (hand and engine driven) such as:

- Endodontic instruments (files, broaches, reamers, enlargers, endodontic burs, ultrasonic inserts);
- Rotary cutting instruments (Diamond burs, tungsten carbide burs, stainless steel drills, carbon steel burs).

Root canal filling instruments (Pluggers, spreaders, compactors).

Supports, kits and instrument organisers.

Hand instruments and clamps.

A2. Implantable devices

Dentinal and radicular posts made of steel, titanium and glass fibers.

Supports, kits and organiser systems for posts.

A3. Contra angle

B. Filling material: Only chemical disinfection (no sterilisation)

Gutta percha, Thermafil obturation devices.

Exclusion

- Equipment such as Motors, Apex locators and other devices with reprocessing procedures included in the individual Direction for Use.
- MTA, Glyde, TopSeal.

General recommendation

- Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer. For all metal instruments, it is recommended to use anticorrosion disinfecting and cleaning agents.
- For your own safety, please wear personal protective equipment (gloves, glasses, mask).
- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterility.
- Limitations and restrictions on reprocessing:
The individual DFU indicates if the useful life of a device might be reduced by the number of reprocessing cycles. Furthermore, the appearance of defects such as cracks, deformations (bent, twisted), corrosion, loss of colour coding or marking, are indications that the devices are not able to fulfil the intended use with the required safety level.
- Single use marked instruments are not approved for re-use.
- The water quality has to be convenient to the local regulations especially for the last rinsing step or with a washer-disinfector.
- Tungsten carbide burs, plastic supports, hand instruments and NiTi instruments are degraded by Hydrogen Peroxide (H₂O₂) solution.
- NiTi Instruments are degraded if immersed more than 5 minutes in a solution of NaOCl at more than 5%.
- Instruments made of aluminium are degraded in presence of caustic soda solutions with mercury salt. Do not use acid (pH < 6) or alkaline (pH > 8) solutions.
- The washer-disinfector is not recommended for instruments made of aluminium, tungsten carbide or carbon steel.

Step-by-step procedure

A. Devices

			A3. Contra Angle				
			A2. Implantable devices				
			A1. Instruments	Following uses			
				First use			
	Operation	Operating mode	Warning				
1.	Disassembling	- Disassemble the device, if required	Silicone stops have to be removed.		X		
2.	Pre-Disinfection	- Soak all instruments immediately after use in a detergent and disinfecting solution combined with proteolytic enzyme if possible.	<ul style="list-style-type: none"> - Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects on instruments). - The disinfecting solution should be aldehyde free (to avoid blood impurities fixation) and without di- or triethanolamines as corrosion inhibitor. - Do not use disinfecting solutions containing Phenol or any products which are not compatible with the instruments (See general recommendations). - For visible impurities observed on instruments a pre-cleaning is recommended by brushing them manually with soft material. 		X		
3.	Rinsing	- Abundant rinsing (at least 1 min)	<ul style="list-style-type: none"> - Use quality water in accordance with local regulations. - If a pre-disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the cleaning. 		X		
4a.	Automated Cleaning with washer-disinfector	<ul style="list-style-type: none"> - Place the devices in a kit, support or container to avoid any contact between instruments or posts. - Put them in the washer-disinfector (Ao value > 3000 or, at least 5 min at 90 °C). 	<ul style="list-style-type: none"> - Discard any instruments with large obvious defects (broken, bent). - Avoid any contact between instruments or posts when placing in the washer-disinfector use kits, supports or container. - Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). - Use only approved washer-disinfector according to EN ISO 15883, maintain and calibrate it regularly. 	X	X	X	X
OR							
4b.	Manual Cleaning or assisted by an ultrasonic device	<ul style="list-style-type: none"> - Place the devices in a kit, support or container to avoid any contact between instruments. - Immerse in the disinfecting solution with cleaning properties, assisted by an ultrasonic device if suitable. 	<ul style="list-style-type: none"> - No visible impurities should be observed on the instruments. - Discard any instruments with large obvious defects (broken, bent, and twisted). - Follow instructions and observe concentrations and time given by the manufacturer (see also general recommendations). - The disinfecting solution should be aldehyde free and without di- or triethanolamines as corrosion inhibitor. 	X	X	X	
5.	Rinsing	- Abundant rinsing (at least 1 min)	<ul style="list-style-type: none"> - Use quality water in accordance with local regulations. - If a disinfecting solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the autoclaving. - Dry on a single use non-woven cloth, or with a drying machine or filtered compressed air. 	X	X	X	
6.	Inspection	<ul style="list-style-type: none"> - Inspect devices and sort out those with defects. - Assemble the devices (stops) 	<ul style="list-style-type: none"> - Dirty instruments must be cleaned and disinfected again. - Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of colour coding or marking) affecting the resistance, the safety or the performance of the instrument or posts. - Protect carbon steel bur with corrosion inhibitor before packaging. - For Contra Angle: lubricate the device with an adequate spray before packaging 	X	X	X	X
7.	Packaging	<ul style="list-style-type: none"> - Place the devices in a kit, support or container to avoid any contact between instruments or posts - Pack the devices in "Sterilisation pouches". 	<ul style="list-style-type: none"> - Avoid any contact between instruments or posts during sterilization. Use kits, supports or containers. - Check the validity period of the pouch given by the manufacturer to determine the shelf life. - Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607. 	X	X	X	X
8.	Sterilization	- Steam sterilisation at: 134 °C / 273°F during 18 min.	<ul style="list-style-type: none"> - The instruments, posts and the plastic supports must be sterilized according to the packaging labelling. - Use only autoclaves that are matching the requirements of EN 13060, EN 285. - Use a validated sterilisation procedure according ISO 17665 - Respect the maintenance procedure of the autoclave device given by the manufacturer. - Use only this recommended sterilization procedure. - Control the efficiency (packaging integrity, no humidity, colour change of sterilisation indicators, physico-chemical integrators, digital records of cycles parameters). - Traceability of procedure records 	X	X	X	X
9.	Storage	- Keep devices in sterilization packaging in a dry and clean environment	<ul style="list-style-type: none"> - Sterility cannot be guaranteed if packaging is open, damaged or wet. - Check the packaging and the medical devices before using them (packaging integrity, no humidity and validity period). 	X	X	X	X

B. Filling material			
	Operation	Operating mode	Warning
1.	Disinfection	- Immerse the obturation devices in NaOCl (2,5% at least) during 5 mn at ambient temperature.	- Do not use disinfecting solutions containing Phenol or any products which are not compatible with the treated filling material. (See general recommendation).

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